# Chemical range 2/2

## Instructions for use

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## GB / CHEMICAL RANGE FIELD OF APPLICATION

The CE marking on these products means that they meet the requirements of EU Regulation 2016/425 on Personal Protective Equipment concerning protection, comfort and strength. Gloves meet the requirements (mechanical, comfort, robustness and protection against the risks claimed) of the PPE regulation 2016/425. The CE marking is issued by CTC (reference number 0075). For further information on the CE marking, please refer to the PPE regulation (EU) 2016/425 on personal protective equipment, as amended to apply in GB, the UKCA marking was issued by SATRA technology centre Ltd (AB0321). • Gloves for protection against chemicals such as acids, bases, detergents, alcohols, ketonic solvents, petroleum solvents, aromatic and chlorinated solvents within the limits of the restrictions specified in the chemical resistance table and/or against microorganisms and/or to provide thermal protection (hot or cold) and/or against radioactive contamination and/or against mechanical risks. • Resistance to the penetration was evaluated under laboratory conditions and only concerns the test specimen under investigation. The gloves do not contain substances at levels such as are known or suspected to have harmful effects on the hygiene or health of the user under foreseeable conditions of use. The 651 glove complies with the EN16350:2014 standard. Maximum recorded RV : Palm 2, 15 x 10<sup>5</sup> Q / Back 1, 18 x 10<sup>5</sup> Q / Cuff 1, 18 x 10<sup>5</sup> Q. • The person wearing the electrostatic dissipative protective gloves must be suitably connected to the ground, for example by wearing appropriate footwear. The electrostatic dissipative protective gloves must not be removed from their packaging, opened, adjusted or removed in flammable or explosive atmospheres, or when handling flammable or explosive substances; The electrostatic properties of the protective gloves may be adversely affected by ageing, wear, contamination and degradation. They may not be suitable for use in oxygen-enriched flammable atmospheres where further assessment. • Gloves with reference 493 for handling phytosanitary products meet the minimum level 2 permeation requirements for the following products: - Isopropanol (Permeation = 6, Degradation = -13) - Cyclohexanone (Permeation = 3, Degradation = 63) - Xylène (Permeation = 2, Degradation = 54) • The permeation levels obtained do not reflect the actual duration of protection in the workplace, nor the differentiation between mixtures and pure chemicals. • The chemical resistance was evaluated under laboratory conditions from samples taken only from the palm (except where the length of the sleeve of the glove was greater than or equal to 400 mm as also checked) and only concerns the chemical subject of the test. It can be different if it is used in a mixture. • Gloves giving protection from radioactive contamination do not protect from ionising radiation and use of the gloves does not prevent the user from being irradiated. • The gloves are not designed to be used in containment enclosures. They may be used as an under-glove for handling waste and for routine cleaning work. • Gloves containing natural latex: avoid contact with oils and petroleum, aromatic or chlorinated solvents. • For nitrile gloves: avoid contact with ketones and organic nitrogen products. • For neoprene gloves: avoid contact with aromatic and chlorinated solvents. • For PVC gloves: avoid contact with ketones and aromatic or chlorinated solvents. • For Butyl gloves: avoid prolonged contact with aromatic solvents and hydrocarbons. • For fluoroelastomer gloves: avoid contact with ketones and acetates. • For category III gloves - Protection against fatal or irreversible hazards: Module D, monitored by ASQUAL - 0334. Module C2, monitored by CTC - NB 0075. • For category III gloves - Protection against fatal or irreversible hazards: Module D/C2, monitored by SATRA UK - 0321

## INSTRUCTIONS FOR STORAGE AND USE

It is recommended that you pre-test the gloves as the actual workplace conditions of use may differ from those of the CE type tests (in particular mechanical and/or chemical), according to temperature, abrasion and degradation. • When used, protective gloves may offer less resistance to dangerous chemicals due to the alteration of their physical properties. The movements, rips, friction or degradation caused by contact with chemicals, etc. can significantly reduce the actual useful life. • For corrosive chemicals, degradation may be the most important factor to be considered when choosing chemical resistant gloves. Before use, it is recommended to inspect the gloves to ensure they do not show any defect or imperfection. • Store the gloves in their original packaging away from light, heat and humidity, in particular, neoprene gloves should be stored at a temperature above 5°C. The integrity of the gloves shall be checked before use (presence of holes, cracks, tears, etc.) and discard any gloves with defects before use. • The design and performance are not significantly affected by ageing if the gloves are stored in the appropriate conditions (humidity, temperature, cleanliness, ventilation, lighting). • Gloves should not be used near machinery due to risk of entrapment. • Thermal protection gloves are designed for limited handling of hot parts up to temperatures of 100°C for a level 1 and 250°C for a level 2. • Do not put the gloves in direct contact with an open flame. The thermal performance level only applies to the coated parts of the glove. • Nitrile or latex coated gloves are not recommended for use by those sensitive to thiocarbamates and/or thiazoles. • For gloves coated in natural latex or natural blended latex; not recommended for use by those sensitive to natural latex and thiuram. • Make sure that hands are clean and dry before putting the gloves on. • Clean the gloves before removing them. These gloves are not machine washable. Use with solvents (diluents, etc.); wipe with a dry cloth. Use with detergents, acids or alkaline products: rinse thoroughly with running water and wipe with a dry cloth. Use with paints and inks: clean with a cloth soaked in a suitable solvent, then wipe using a dry cloth. Use with phytosanitary products: Using water, immediately wash gloves contaminated with undiluted product. Then add the resulting rinse water to the container holding the chemical spray. • Caution: improper use of the gloves or cleaning them in a way that is not specifically recommended can alter their performance levels. • Ensure the inside of the gloves is dry and that they are in good condition before reusing them. • For more information about the performance levels, chemical resistance and usage of the gloves, please contact your distributor or MAPA PROFESSIONAL Technical Customer Support. • Information leaflet and EU/UKCA Declaration of Conformity can be downloaded from [www.mapa-pro.fr](http://www.mapa-pro.fr)

MAPA SPONTEX UK Ltd, Berkeley Business Park Wainwright Road, Worcester WR4 9ZS  
T: (44) 1 905 450300 - F: (44) 1 905 450350 - DG 1 905 450360. [www.mapa-pro.uk](http://www.mapa-pro.uk)

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MAPA S.A.S. Défense Ouest - 420, rue d'Estienne d'Orves F - 92705 COLOMBES Cedex
MAPA GmbH Industriestraße 21-25 D - 27404 Zeven
Mapa Spontex Ibérica S.A.U. Llacuna, 161 - Planta 3ª, Módulo D - 08018 BARCELONA