- DE Benutzerinformation
- EN User information
- CZ Informace pro uživatele
- DK Brugervejledning
- FR Notice d'utilisation
- GR Πληροφορίες για τον χρήστη
- NL Gebruikersinformatie
- PL Informacja dla użytkownika

- DE Untersuchungshandschuhe
- EN Examination gloves
- CZ Výměňovací rukavice
- DK Undersøgelseshandsker
- FR Gants d'examen
- GR Γάντια Εξέταστικά
- NL Onderzoekshandschoenen
- PL Rękawice diagnostyczne

- DE Einmalschutzhandschuhe
- EN Disposable protective gloves
- CZ Jednorázové ochranné rukavice
- DK Engangs- beskyttelseshandsker
- FR Gants de protection à usage unique
- GR Προστατευτικό γάντι μιας χρήσης
- NL Wegwerp beschermhandschoenen
- PL Jednorazowe rękawice ochronne

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**GB User Information**  
Protective Gloves against irreversible Risks / Category III

ISO 374-1 Type A ISO 374-1 Type B ISO 374-1 Type C

TYPE A Permeation min. 6 chemicals, level 2 (table 1)  
TYPE B Permeation min. 3 chemicals, level 2 (table 1)  
TYPE C Permeation min. 1 chemicals, level 1 (table 1)

ISO 374-5:2016  
Protective gloves protecting from viruses, bacteria and fungi

ISO 374-5:2016  
Protective gloves protecting from bacteria and fungi

If not checked for protection against viruses the following warning must be added: "Not checked against viruses".

| Table 1                  | Time     | Protection index |
|--------------------------|----------|------------------|
| A methanol               | >10 min  | 1                |
| B acetone                | >30 min  | 2                |
| C acetonitrile           | >60 min  | 3                |
| D dichloromethane        | >120 min | 4                |
| E carbon disulphide      | >240 min | 5                |
| F toluene                | >480 min | 6                |
| G dimethylamine          |          |                  |
| H tetrahydrofuran        |          |                  |
| I ethyl acetate          |          |                  |
| J n-heptane              |          |                  |
| K sodium hydroxide 40%   |          |                  |
| L sulphuric acid 96%     |          |                  |
| M nitric acid 65%        |          |                  |
| N acetic acid 99%        |          |                  |
| O ammonium hydroxide 25% |          |                  |
| P hydrogen peroxide 30%  |          |                  |
| S hydrofluoric acid 40%  |          |                  |
| T formaldehyde 37%       |          |                  |

**EN 374-4:2013**

The degradation (in %) indicates the change in puncture resistance of the gloves after exposure to the respective challenge chemical.

**Compliant with PPE Regulation (EU) 2016/425, EN 420 and ISO 374.**  
The gloves are classified as Personal Protective Equipment (PPE) Category III according to PPE Regulation (EU) 2016/425 and have been shown to comply with this regulation through the applicable harmonised European standards.

XXXX = Identification number of notified Body responsible for the EU type examination and supervising ongoing conformity. (see box)

**Storage/expiry date:** Store gloves in the original packaging on a cool and dry place. Never place heavy items on top of it. Keep away from direct sunlight. The storage place shall not contain any equipment that generates ozone (laser printers, copiers etc.) or others that give rise to electrical sparks or discharges.

The actual date of expiry cannot be stated as it depends on the general conditions when being used. An individual risk assessment is required in any case.

Expiry date when stored as intended.

**Use and control:** Protective gloves should only be used for the intended application and in the correct size. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. Before use the gloves should be checked for any defects such as holes, cracks or further faults like discolourations. Defective gloves may not be used under any circumstances. The breakthrough times and permeation levels are based on laboratory measurements and are determined using samples from the palm of the hand. The actual duration of protection of a glove with a certain substance can differ significantly due to the conditions of use (temperature, abrasion, elongation). For aggressive chemicals, degradation (change in mechanical properties) can be the most important factor to consider in selection of chemical resistant gloves.

**Cleaning and disinfection:** There is no cleaning or disinfecting provided for these gloves. If carried out anyway, this will be at your own risk.

**Disposal:** Used gloves must be disposed of after contact with chemicals according to the disposal regulations of the chemical. Unused gloves may be disposed of together with household waste.

**Ingredients / Hazardous ingredients:**

Some gloves might contain ingredients which are known to possibly cause skin irritations or allergic reactions with sensitised persons. Check warning information on specific packaging carefully.

**Warnings:** Protective gloves are intended for single use only. Do not wear any gloves near rotating machine parts (saw blades, drills etc.). There is the risk of getting caught!

**Declaration of conformity:** www.igefa.de

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